


Phonak Marvel™ Receiver and Coupling Guide

Receivers


Standard (S)
46/111 (2cc)



Mild to moderately-severe hearing loss


Preferred Marvel receiver

Medium (M)
50/114 (2cc)




Mild to severe hearing loss

Power (P)
58/122 (2cc)



Mild to severe hearing loss

UltraPower (UP)
67/130 (2cc)
Compatible with cShell only



Moderate to profound hearing loss

Phonak Audéo™ M



312



R



312T



13T



RT

Available Fall 2019

Coupling Options



Cap Dome



Open Dome



Vented Dome



Power Dome



CeruShield™



SlimTip
(Acrylic)



RIC SecureFit
(Silicone)



SlimTip
Titanium



cShell